

Title: Multicenter Systems Analysis of Human Blood Reveals Immature Neutrophils in Males and During Pregnancy

Authors: Blazkova et al

Publication: Journal of Immunology, 2017

Data Source: Supplementary Table 3

Antibody panel used in mass cytometry analysis of samples banked Smart Tube Proteomic Stabilizer PROT1

| Antigen | Symbol | Atomic Mass |
|--------------------|--------|-------------|
| CD235ab | In | 113 |
| CD61 | In | 113 |
| CD45 | In | 115 |
| CD66 | La | 139 |
| CD7 | Pr | 141 |
| CD19 | Nd | 142 |
| CD45RA | Nd | 143 |
| CD11b | Nd | 144 |
| CD31 | Nd | 145 |
| CD8a | Nd | 146 |
| CD123 | Nd | 148 |
| TCR $\gamma\delta$ | Sm | 152 |
| CD10 | Gd | 156 |
| CD161 | Gd | 157 |
| CD33 | Gd | 158 |
| CD116 | Dy | 162 |
| CD16 | Ho | 165 |
| CCR9 | Er | 168 |
| CD25 | Tm | 169 |
| CD3 | Er | 170 |
| CCR7 | Yb | 171 |
| CD15 | Yb | 172 |
| CCR2 | Yb | 173 |
| HLA-DR | Yb | 174 |
| CD14 | Yb | 175 |
| CD56 | Yb | 176 |