

Title: Mapping the Fetomaternal Peripheral Immune System at Term Pregnancy
 Authors: Fragiadakis et al
 Publication: Journal of Immunology, 2016
 Data Source: Supplemental Table 1

Antibody panel used in mass cytometry analysis of samples banked Smart Tube
 Proteomic Stabilizer PROT1

Antigen	Symbol	Atomic Mass	Clone	Supplier
CD235ab	In	113	HIR2	Biolegend
CD61	In	113	VI-PL2	BD
CD45	In	115	HI30	Biolegend
CD66	La	139	CD66 α -B1.1	BD
CD7	Pr	141	M-T701	BD
CD19	Nd	142	HIB19	Fluidigm
CD45RA	Nd	143	HI100	Fluidigm
CD11b	Nd	144	ICRF44	Fluidigm
CD4	Nd	145	RPA-T4	Fluidigm
CD8a	Nd	146	RPA-T8	Fluidigm
CD11c	Sm	147	Bu15	Fluidigm
CD123	Nd	148	6H6	Biolegend
p-CREB	Sm	149	87G3	CST
p-STAT5	Nd	150	47	Fluidigm
p-P38 (pT184/pT182)	Eu	151	36	CST
TCR γ δ	Sm	152	GL3	Fluidigm
p-STAT1	Eu	153	58D6	Fluidigm
p-STAT3 (pY705)	Sm	154	4/P	BD
p-rpS6 (pS235/236)	Gd	155	N7-548	CST
CXCR3	Gd	156	G025H7	Fluidigm
CD161	Gd	157	HP-3G10	Biolegend
CD33	Gd	158	WM53	Fluidigm
p-MAPKAPK2	Tb	159	27B7	Fluidigm
Tbet	Gd	160	4B10	Fluidigm
FoxP3	Dy	162	PCH101	Fluidigm
I κ B	Dy	164	L35A5	Fluidigm
CD16	Ho	165	3G8	Fluidigm
p-NF κ B	Er	166	K10-895.12.50	Fluidigm
GATA3	Er	167	TWAJ	Fluidigm
CCR9	Er	168	L053E8	Fluidigm
CD25	Tm	169	2A3	Fluidigm
CD3	Er	170	UCHT1	Fluidigm
p-ERK1/2	Yb	171	D13.14.4E	Fluidigm
IgM	Yb	172	MHM-88	Fluidigm
CD127	Yb	173	HIL-7R-M21	BD
HLA-DR	Yb	174	L243	Fluidigm
CD14	Yb	175	M52E	Fluidigm
CD56	Yb	176	NCAM16.2	Fluidigm